

Make 'n Take Sheet

Discover How Easy It Is To Make Your Own Assistive Technology!

ABLE
Tech

Button Pusher/Page Turner

A Button Pusher/Page Turner can provide an easy grip for someone with limited finger and hand control.

Materials:

- | | |
|---|---------------------------|
| 1 - 3/4" to 1" dowel purchased in a hobby, art, or hardware store | Drill with 1/4" drill bit |
| 1 - unsharpened pencil with eraser | Hand saw |

Directions:

1. Measure and saw the dowel in 5 1/2" sections.
2. Drill a 1/4" hole in each section of dowel located at the midpoint and about 1/2" deep.
3. Push the lead end of the pencil into the hold and glue if necessary.
4. Repeat the process for each remaining dowel section.



Oklahoma ABLE Tech, the statewide Assistive Technology Act Program, is funded through the Administration for Community Living of the U.S. Department of Health and Human Services, and maintains coordination and collaboration efforts with partners throughout the State of Oklahoma.

To contact Oklahoma ABLE Tech call 1-800-257-1705