

Make 'n Take Sheet

Discover How Easy It Is To Make Your Own Assistive Technology!

Key Holder

People with impaired hand function achieve a better grip when using this key holder.

Materials:

1 – Wooden dowel or PVC pipe
(approx. 6" length and 1/2"
diameter)

Screwdriver

Key

1 – Screw (1/2")

Directions:

1. Pre-drill a small hole in the middle of the dowel.
2. Line up the middle hole of the key with the pre-drilled hole.
3. Use screwdriver to tighten screw into the dowel or PVC pipe.



Purchase Price: \$5.99



Make-N-Take Price: \$1.25

Helpful Hints: No drill? Ask at a home improvement store to help pre-drill the holes or use a hammer and nail to start a small hole for the screws.

Oklahoma ABLE Tech, the statewide Assistive Technology Act Program, is funded through the Administration for Community Living of the U.S. Department of Health and Human Services, and maintains coordination and collaboration efforts with partners throughout the State of Oklahoma.

To contact Oklahoma ABLE Tech call 1-800-257-1705