


Make 'n Take Sheet

Discover How Easy It Is To Make Your Own Assistive Technology!

ABLE
Tech

Tactile Ruler



Materials:

1 – 12" ruler

12 - small foam circles (6 each of 2 contrasting colors)

Glue

Directions:

1. Use a hole punch to punch holes in craft foam.
2. Glue foam circles to ruler at each number. Alternate colors.

Oklahoma ABLE Tech, the statewide Assistive Technology Act Program, is funded through the Administration for Community Living of the U.S. Department of Health and Human Services, and maintains coordination and collaboration efforts with partners throughout the State of Oklahoma.

To contact Oklahoma ABLE Tech call 1-800-257-1705.